

# SWINGO-T & SWINGO-P

- FR** NOTICE D'UTILISATION - Ancillaire SWINGO-T & SWINGO-P
- EN** INSTRUCTIONS FOR USE - SWINGO-T & SWINGO-P Instruments
- ES** MANUAL DE INSTRUCCIONES - Instrumentos SWINGO-T & SWINGO-P
- DE** BEDIENUNGSANLEITUNG - Zubehör SWINGO-T & SWINGO-P
- GR** ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Βοηθητικό υλικό SWINGO-T & SWINGO-P
- IT** ISTRUZIONI PER L'USO – Strumentario SWINGO-T & SWINGO-P
- PT** MANUAL DE UTILIZAÇÃO - Equipamento auxiliar SWINGO-T & SWINGO-P
- RO** INSTRUCȚIUNI DE UTILIZARE - Material auxiliar SWINGO-T & SWINGO-P

# Ancillaire SWINGO-T & SWINGO-P

## Généralités

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
# SWINGO-T & SWINGO-P Instruments

**Generalities**

OSD recommends that all medical staff in contact with instrument take knowledge of these instructions for use. The using of the material requires anatomy knowledge, biomechanical and spine surgery knowledge, material can be used only by a qualified surgeon who practices this surgery taking account science and surgical art data available at this moment.

**Description**

Instrument is composed of the instrumentation (instrument material), storage and a container. It is necessary to take knowledge of the surgical technique and in case of any doubt; contact OSD Company agent's to obtain all information necessary for using of the instrument and associated implants. Besides, some instruments are not delivered with standard set, they can be provided on request. In no way, an instrument or an instrument component can be implanted. Instrument material can be made in stainless steel, in titanium alloy, in polymer and in silicon. The materials are biocompatible, besides all modification provided to devices may call into question their technics and biologics.

**Indications**

OSD material instrument has been developed especially to facilitate the implantation or the ablation of SWINGO implants and must be used exclusively for this use. The health institute is responsible to proceed to the pre-disinfection, cleaning and the sterilization of the instruments after its using according to the validated methods. The advices attached not are a substitute for current sanitary rules: norms, good practice, handbook, national advices, ministerial texts, etc.

**Contraindications**

Allergy known to one of the material component (in particular to nickel). The contraindications are associated mostly for SWINGO (instructions for use SWINGO IMPLANTS : Rev\_20190501\_SWGS\_IFU) OSD accepts no responsibility in case of indications and contraindications non-respect concerning SWINGO instrument and implants.

**Functionality inspection and verification of SWINGO instruments**

As a precautionary measure against possible patients and users injuries, a control must be made every time at the reception of the instrument and before each using. All breakdowns can cause serious risks. At the reception, instruments set must be controlled with precision. The reference of the instrument is written on the device. Furthermore, it is important to verify the availability of all devices, if need be different sizes constituting set instrument, from instructions and expressed needs by the surgeon. A list of delivered material is attached to the instrument shipment.

When an instrument is associated to an other during the surgical process, it is necessary to verify that the device is correctly assembled. Verify the action of moving part to insure good function, precisions and using modalities are available on OSD website: <http://osdevelopment.fr>

Inspect carefully each device to insure that all visible contamination has been eliminated or all damage and/or excessive wear are identified.

In case of visible contamination, it is obligatory to renew cleaning and disinfection process.

All damaged device must be isolated and returned at OSD.

**Side effects**

Allergy to one of the material component.

Side effects SWINGO implants (instruction for use SWINGO IMPLANTS : Rev\_20190501\_SWGS\_IFU).

**• Warnings and precautions**

Precautions must be taken during the handling of devices in particular in case of sharpened or edged devices.

In case of an instrument is delivered alone. At reception it is necessary to remove plastic bag before cleaning and sterilization.

The articulated instruments must be opened, to improve the washing efficacy.

In order to ensure the functionality and the security of the material, the instructions attached must be respected:

- Opening of the packaging must be realized near surgical site.
- Medical devices must be handled with care.
- The using of brushes, metal stamp scrub or all tools which may alter the instrument during manual cleaning process is forbidden.
- Favour the using of soft nylon brushes.
- The using of chemical products with chlorine, acids, and organics solvent or ammoniacal is forbidden.
- The using of product with soda is forbidden for the container.
- The using of acid for the neutralization of alcalin residue at the end of automatic washing cycle on aluminum trays, and instruments comportsing polymer elements is forbidden.
- In general, it is necessary to pay attention to the different elements in case where a part was lacking, lost or broken; we invite you to contact OSD.
- Not air dry the instruments after using, the soiled devices which have been dried are harder to clean.
- Not use mineral oil or lubricant with silicon because they covered micro-organisms, prevent the direct contact of the superficies with the vapor and are hard to eliminate.
- Not bring modification to the instruments.
- The saline solution and cleaning agents/disinfection containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and must be not used.
- Cleaning agents and enzymatic with neutral pH are recommended.

**Note: Orthopedic acts are not considered like risk at the ATNC. A total inactivation with molar soda (1N) or sodium hypochlorite with a concentration of 2% of active chlorine not is reserved to the instruments used previously on suspect or suffering EST before invasive act.**

**• Precaution / retreatment**

If need be, it is obligatory to remove the surplus of organics liquids and tissues of the instruments with a lint-free disposable wipe.

It is recommended to not air dry the saline solution, blood, osseous fragments, tissue or other organic debris on the instruments before cleaning them.

Place the instruments in distilled water bath or in a tray covers with damp towels.

Equipment for personal protection must be worn during the handling of the contaminated devices.

**• Pre-Disinfection**

The pre-disinfection has to aim to reduce the micro-organism population and to facilitate the ulterior cleaning. It has also to aim to protect the staff during the handling of the instruments and avoid the contamination of the environment. All re-usable device shall undergo an immediately pre-disinfection or be treated without delay in disinfect scrubbed.

WARNING: Containers and storage trays must not be in prolonged contact with this type of decontaminant; clean the soiled zone and rinse immediately.

- **ATNC**

It is necessary to respect the dispositions in:

- French instruction DGS/RI3/2001/449 of 1<sup>st</sup> December 2011 pertaining to the updating advices aiming to reduce disease risks of non-conventional transmissible agents during invasive act.
- French circular DGS/SD5C/DHOS/2005/435 of 23th September 2005 pertaining to the advices for medical devices treatment used on subjects who have received labile blood products (PSL) coming from donors retrospectively affected on cases of variant Creutzfeld-Jakob (vMCJ) disease.

## **Cleaning, disinfection and sterilization instructions:**

- **Cleaning / disinfection**

Cleaning and disinfection are obligatory before sterilization of re-usable material according to the indicated process. This step has to aim to remove all soiling and residue on the instruments cause to the surgery, in order to reduce the number of viable micro-organism on the instruments before sterilization.

Cleaning and disinfection methods common used are: manual cleaning and automatic cleaning.

The choice of implementation is under responsibility of the health institute.

In order to ensure the cleaning and disinfection efficacy, please respect the instructions for use of the cleaning product.

During this step, a special attention should be given to sensible zone and/or which has a difficult access so that the product and the cleaning action are optimal.

In first time, it is important to remove the liquids excess and organics tissues with a lint-free disposable wipe. And place the devices on a distilled water tray or cover with a damp towel.

### *Cleaning/automatic disinfection process:*

Place the instruments in the washer/disinfector adapted basket, according to the instructions of the manufacturer and treat via the recommended cycle for OSD instruments and described below:

- **Pre-cleaning**

- Disassembly of the instruments, if necessary.
- Immersion during 5 min at least on cold tap water.
- Brush the superficies with a nylon brush until the disappearance of visible residue.
- Rinsing during 20 seconds under pressure jet (between 3 and 4 bar or 43,5 and 58 psi), insisting on cannula.

- **Washer/disinfector**

Use a washer/disinfector type MIELE G7735 CD as well as its program Vario TD, according to these steps:

- 2 minutes of pre-cleaning on cold water.
- Flush.
- 5 minutes of cleaning on demineralized water at 55°C and 0,5% of soft alcalin detergent (for example : NeodisherMediClean - Dr.Weigert).
- Flush.
- 3 minutes of rinsing and neutralization with demineralized water.
- Flush.
- 2 minutes of final rinsing with demineralized water
- Flush.
- Drying, dry cloth, lint free, clean and sterile.

Verify that the instruments and the lumens have not residue.

If necessary, repeat the cleaning cycle.

After cleaning/disinfection, the disassembled instruments must be reassembled and inspected visually.

- **Inspection before sterilization**

Inspect with attention each device to ensure the elimination of all visible contamination. If necessary, repeat the cleaning and disinfection process.

Furthermore, verify the absence of the instruments deformation.

- **Sterilization**

The obtention and maintain of the sterile state of instruments are under the responsibility of the health institute. The sterilization with vapor/moist heat is the recommended method.

Furthermore, OSD commercial service is here to help you for all complementary information.

The dispositions to take for the protection of sharped crest or potentially injuring instruments must be taken in consideration by the health institute.

OSD recommends respecting the recommendations of the manufacturer concerning the sterilizer, which must be every time followed. During the sterilization several instruments sets on the sterilization cycle, it is important to ensure do not exceed the maximal load indicated by the manufacturer.

Instruments sets must be prepared and packaged correctly in trays and/or boxes to permit to the vapor to penetrate and enter in direct contact with all the superficies.

Refer to the below information for the sterilization parameters of validated cycles by OSD in order to give an insurance level of sterility (NAS) of  $10^{-6}$ .

Sterilization methods to ethylene oxide or to plasma not must be used.

Sterilization parameters to the validated vapor:

Type of cycle: air removal - Temperature 134° C (273.2°F)

Exposition time: 18 minutes – Drying time: 20 minutes.

## **REPARATIONS**

All device damaged must not be used and must be returned for reparation at OSD.

In order to protect the personal responsible of the reparation, instruments can not be returned for reparation only if the Decontamination/Cleaning/Sterilization processes are respected. If the instruments are soiled or if modifications are noted on the products to repair, OSD reserves the right to repair these medical devices to your cost or to refuse the entire reparation.

OSD is not responsible in case of modification or reparation effectuated on the devices by the health institute.

## **ELIMINATION**

If need be, the device can be eliminated according to the good practice of infectious elimination waste. There is no special risk or unusual associated to the device elimination.

## **SHIPMENT RETURN**

Health institute must disinfect, clean, sterilize and package correctly the instruments (complete box or isolated instrument) before returning them at OSD.

Furthermore, when the health institute has not the medical device property, it accepts invoicing and takes all the necessities for the payment in the following cases:

- When the medical device is destined to be used in circumstances where there is suspicion of ATNC contact (patient or surgical act with risk);
- When the destruction of the medical device is imposed as part of the observation of sanitary rules and also if OSD instructions are not respected, instructions mentioned on this instructions for use concerning the storage, decontamination, cleaning and sterilization and inducing a material deterioration.

**GUARANTEE**

OSD Company guarantees that its products delivered are compliant to the associated specifications. Furthermore, this guarantee can not be applied in case of bad using, or bad handling of SWINGO device. OSD is not responsible concerning the additional using accessories used during the implantation of SWINGO device. All defective instruments must be returned at OSD to be replaced or destroyed. This instructions for use must be accessible and kept close to the instrumentation, OSD should not be responsible of damages resulting a bad using and/or device handling and in case where repair and/or modification should be given to the instrument by a non-authorized person.


OSD Company does not guarantee the clinical result, knowing that residual risks can compromise the success of the surgery.

To ensure an optimal using of devices, it should be respected the following symbols:

Symbols used						
						
Consult instructions for use	Manufacture	Not sterile	Do not use if package is damaged	Catalog Number	Keep away from sunlight	Keep dry

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